sartorius stedim
biotech

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Vers. 11 | 2011

Table of Contents

1	Use of the Calibration Unit	3
2	Description of the Calibration Unit	3
3	Specifications	4
4	Equipment Supplied	5
5	Operation of the Calibration Unit	5
6	Calibration of the MD8 airscan	7
7	Calibration of the AirPort MD8	8
8	Care and Maintenance	9
9	Information and Instructions on Disposal	9

sartorius

Declaration of Conformity CE
According to Council Directives 89/336/EEC and 73/23/EEC
(revised by Directive 93/68/EEC)

Company	Sartorius AG
Address	Weender Landstr 94-108 D-37075 Göttingen Germany
	<p>We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius AG.</p>
Designation of Device	Calibration Unit for MD8
Cat.-No	16756
Relevant directives of the EU	EU Directive on Electromagnetic Compatibility (89/336/EEC) in Version 92/31/EEC EU Low Voltage Equipment Directive (73/23/EEC) revised by Directive 93/68/EEC
Applied harmonized standards	<p>EN 61326-1 : Electrical equipment for measurement, control and laboratory use. Limitation of emissions: Residential areas B Defined immunity to interference: Industrial areas, continuous unmonitored operation</p> <p>EN 61010 : Safety requirements for electrical equipment for measurement, control and laboratory use</p>
	Test protocol No. 289, Sartorius AG, Abt. EMC
Date and Signature	 
Function of the Signatory	<p>V. Niebel Dr. Susanne Gerthausen</p> <p>Senior Vice President Quality Management Representative</p>

1 Use of the Calibration Unit

The Sartorius 16756 calibration unit is a portable flow rate measuring device that was specially developed for flow rate calibration of the MD 8 airscan and AirPort MD8 air samplers.

2 Description of the Calibration Unit


Fig. 1

The unit consists of a stainless steel housing with an air inlet port for attachment of a 17655 filter holder or adapter for the gelatin filter units on the front panel (Fig. 1, see a) and an air outlet port on the rear panel (see b) for attaching an 17085 hose for connection to the air sampler (bayonet-design quick connect). Located on the right side (Fig. 1, c) is the connector for the battery charger | power supply unit as well as an opening (d), through which the reset button located inside can be operated, for example, by inserting a ballpoint pen. This reset button is used to switch off the unit by circumventing the automatic shutoff (the calibration unit is not equipped with an OFF switch; the calibration unit switches off automatically after 1 minute).

The start button (e) is located on the upper side of the calibration unit to the right of the LCD display (Fig.1). This button is used to switch the unit on and is then in the ready mode.

The accuracy of the flow rate measurement depends on the horizontal alignment of the calibration unit. For this purpose, a display is located on the right on the upper side (level indicator, f) to allow accurate horizontal alignment and is the reason why the underside of the unit is equipped with 4 adjustable feet (Fig. 1, g).

The RS 232 serial interface (9-pin sub D male connector, h) is located on the lower left of the upper panel, and used for porting peripheral devices, like a PC with a terminal program, for reading and storing data.

2 foldable carrying handles are located on the right and left side panels (Fig.1, i)

3 Specifications

Important Notes

Caution!

- Do not expose the calibration unit to strong vibrations.
- To prevent contamination with airborne particles, do not operate the calibration unit without a filter in place.
- The maximum flow rate is 16 m³/h. This rate may not be exceeded over long time periods.
- Only connect the calibration unit to MD8 airscan | AirPort MD8 air samplers with a regulated air flow rate. Never connect the calibration unit to an air sampler with unknown flow rate or you may possibly overload the flow rate measuring device.
- The unit can be operated independent of mains power supply. When operating the unit on mains power supply however, make sure that voltage rating is correct. If the unit is operated independently of mains power most of the time, it should not be recharged using the power supply unit until the prompt appears in the display. This minimizes a drain on battery capacity caused by the battery's memory effect. When the unit is permanently operated on the mains power supply, it is recommended to leave the unit connected to the mains power supply even when it is turned off because the built-in automatic charging system switches over to a very low charge conservation current during these times.
- Do not expose the unit to extreme fluctuations in temperature or humidity during operation or storage. If fluctuations in temperature and humidity cannot be avoided, allow 2 hours for acclimatization.
- Use other than for the intended purpose, or improper operation, frees the manufacturer from all liability.

Measuring range:	min. 1 m ³ /h	max. 16 m ³ /h
Max. error of the measuring system:	±2% (1 m ³ /h–16 m ³ /h)	
Measuring signal:	2500 pulses = 1 m ³	
Measuring range of the frequency measurement:	0.6 Hz – 12 Hz	
Accuracy:	±0.005 Hz	
Automatic shutoff after 1 minute at a frequency of:	< 0.4 Hz	
RS232 serial interface:	19200 baud, 8 data bits, 1 stop bit, no parity	
Current consumed by measurement electronics:		
Quiescent current:	120 µA	
Operating current:	225 mA	
Battery pack:	NiCd battery pack 9.6 V	
Capacity:	1000 mAh	
Dimensions of battery pack:	96 × 28 × 28 mm	
Weight of the battery pack:	205 g	
Hours of operation with fully charged battery pack:	approx. 4 h	
Time it takes to load a completely discharged battery pack:	approx. 10 h	
Battery monitoring*:		
Automatic charging module	When the calibration unit is supplied by mains power, the battery pack is charged through an automatic charging module. Once the preset end-of-charging voltage (10.6 V) is reached, the charging current (100 mA) is automatically reduced. This prevents the battery pack from overcharging.	
Power supply requirements:	14.5 V DC 530 mA	
Allowable ambient temperature of the calibration unit:	0°C +40°C	
Type of protection:	IP 40 during operation with filter holder, connector and PVC hose	
Numerical display of flow rate on the LCD display (in m ³ /h and l/min)	Alphanumerical display of 4 × 20 pixels	
Weight of the calibration unit without accessories:	Approx. 11 kg	
Dimensions:	Length = 300 mm (without filter holder) Width = 390 mm with handles Height = 182 mm min., 200 mm max. (adjustable feet)	
Type of connector on the inlet and outlet ports for the filter holder and MD8 airscan/AirPort MD8:	Bayonet-design quick connect	

* The calibration unit can be run permanently on the battery charger without overloading the rechargeable battery pack. The unit can also be used when recharging (run on mains power supply).
We recommend leaving the calibration unit connected to the battery charger | power supply unit when it is not being used to prevent the rechargeable batteries from becoming completely discharged.

4 Equipment Supplied

4.1 Equipment Supplied for 16756

- Calibration unit
- Battery charger | power supply unit specific for the country in which it is used
- Filter holder 17655
- 1 set of connectors 17657, consisting of 17658 (for connecting the air sampler to the hose) and 17659 (for connecting the hose to the calibration unit)
- 1 PVC hose (2 m) 17085

Please use either gelatin filter units (order no. 17528-080) or gelatin disc filters (order no. 12602-080).

5 Operation of the Calibration Unit

(Description is based on the example of use with aluminum filter holder 17655).

5.1 Preparation of the Filter Holder

The filter holder is identical to the one on the MD8 airscan | AirPort MD8 air samplers. See the operating instructions for the MD8 airscan | AirPort MD8 air samplers regarding preparation of the filter holder, placing the filter in the filter holder and closing the filter holder.

5.2 Attaching the Filter Holder to the Calibration Unit

Fig. 2

Place the filter holder prepared with the gelatin disc filters (or the aluminum adapter for attaching gelatin filter units) over the air inlet port and attach the filter holder (or the aluminum adapter with gelatin filter unit) to the calibration unit by turning it clockwise (Fig. 2).

5.3. Connecting the Calibration Unit and the Air Sampler

- When connecting the calibration unit to the MD8 airscan | AirPort MD8 consider the direction of air flow!

Fig. 3

Fig. 4

Fig. 5

Turn the 17658 connector (hose | air sampler) clockwise to attach it to the air intake port (Fig. 3) of the air sampler (during calibration, no filter holder is attached to the air sampler). Attach the hose to the end of the connector (Fig. 4) and attach the 17659 connector (hose | calibration unit) to the free end of the hose. Now connect the 17659 connector to the air outlet port on the calibration unit by turning it clockwise (Fig. 5).

Make sure that the connections and seals between calibration unit and air sampler are tight. Figure 6 shows how the final setup looks when it is ready for a calibration to be carried out.


Fig. 6

5.4 Switching on and Operating the Calibration Unit

Press the green start button to switch the calibration unit on. The following message will appear in the display for 5 seconds: „Sartorius MD8 Calibrating System, Version 3.12“ (Fig. 7).


Fig. 7

Note

If the battery voltage is less than 8 Volt, the calibration unit will shut itself off immediately. In such an event, the power supply must be connected and the calibration unit recharged for at least 10 hours.

After 5 seconds, “Sartorius MD8 Calibrating System start air sampler!” will appear in the display (Fig. 8).


Fig. 8

Within 60 seconds of now, the calibration should be started. If the calibration is not started, the unit will switch itself off after 60 seconds. The unit can be switched back on again by pressing the green start button again.

Note

The unit can also be switched off by pressing the reset button. After reset is pressed, the program in the calibration unit will be reset to the default values.

If the calibration unit identifies any air flow during this 60-second waiting time, the air flow rate will be indicated on the display after a delay of 3 seconds:

“Sartorius MD8 Calibrating System”
 m³/h xx
 l/min yy

For reading the calibration data, the display will indicate the rate in m³/h and l/min at the same time (Fig. 9)


Fig. 9

The calibration individual steps are described in Sections 6 (for MD8 airscan) and Section 7 (for the AirPort MD8).

Note

If during operation (calibration procedure) the battery pack voltage drops below 8.8 V, the following message will appear in the display

“Battery is Empty, Plug Power Supply”

The two lines in the display below this will continue to indicate the air flow rate in m³/h and l/min.

In this case, the power supply included in the equipment supplied must be connected to the mains power.

The warning message in the display disappears after a brief period of external power supply and “Sartorius MD8 Calibrating System” appears again.

If the unit is not connected to an external power supply (mains power), the calibration unit will switch itself off automatically (once it reaches its discharge voltage). This prevents the battery pack from discharging too low.

6 Calibration of the MD8 airscan:

(See also Operating Instructions for MD8 airscan Air Sampler)

A complete reference curve for the flow rate sensor is stored in the MD8 airscan software, which can be brought into line with each instrument by the 3-point calibration. In the 3-point mode, the air flow rates at the 3, 5 and 7 m³/h settings are successively run and compared with the values measured by the calibration unit 16756.

Calibration can be interrupted at any time by pressing Stop. When the unit becomes overheated, calibration is interrupted automatically.

The points measured and recorded up to this interruption will be remain stored and the unchanged points are also stored.

You can call up the calibration menu when the MD8 airscan is in the initial starting position, i.e. the message "Sartorius MD8 ready to start" is displayed, by simultaneously pressing the \uparrow and \downarrow keys for 3 seconds.

The calibration function is shown in the display by the message: "- Calibration - 3 point (Start)".

When Start is pressed with the unit in this position, the 3-point calibration mode is started.

6.1 Calibration Procedure

(After connection of the calibration unit is completed as described in Section 5).

- a) The MD8 airscan moves to the first measurement point, and continually displays the nominal and actual values
(The message in the display reads: "Nominal 3.0 m³/h Actual: xx m³/h").
- b) As soon as the nominal value is reached and stabilized, which can be recognized from the actual value display, press the parameter key \rightarrow .
- c) The MD8 airscan now continues to run uncontrolled, i.e. the set point of the motor does not change, so that the air flow rate remains constant. In this mode, the display reads:
"Nominal 3.0 m³/h real: x.xx m³/h".
- d) As soon as the connected calibration unit 16756 has stabilized and displays a constant value, set this value on the MD8 airscan by pressing \uparrow or \downarrow . Then confirm and store it with Enter.
- e) If this is not the last calibration point, the next one is automatically activated and processed as described under a-d.
After the 7 m³/h value has been calibrated, the calibration is finished. This is shown by the message "Calibration completed (Enter)".
- f) Confirm by pressing Enter.
The MD8 airscan returns to the initial starting position and is ready to start.

Note:

If the calibration is interrupted ahead of time, either the message "User break! (press Enter)" or "Overheated device (Enter)" will be displayed. Confirm this message by pressing Enter. The MD8 airscan returns to the initial starting position.

The altered calibration curve is stored immediately when the calibration has been finished or interrupted.

6.2. Resetting of Default Values

If, when calibrating the MD8 airscan, data were entered which are so absurd that they make it impossible for the instrument to be controlled correctly by the phase controller, it may be very difficult to recalibrate correctly.

Even though this case is very unlikely, because the consequences would be serious, a hidden option has been implemented in the MD8 airscan, allowing it to be reset to its default values.

Procedure:

Switch off the MD8 airscan. While holding down both the \rightarrow and Stop keys simultaneously and switch the MD8 airscan on. Do not release the keys until the message "Self-testing" is displayed. Now press the \uparrow key while this message is displayed.

Pressing of the \uparrow key must be made while the "Self-testing" message is displayed as only then will the display change to "Delete data? <Start> = yes". By now pressing Enter you confirm this and all data will be overwritten to the default settings. A self-test is then re-started.

If you press any other key, instead of Enter, the MD8 airscan will enter its standard operating mode without altering the data, in other words, the unit is reset to the default values.

7 Calibration of the AirPort MD8:

7.1 Connecting the AirPort MD8 to Calibration Unit 16756

To run the calibration, connect the calibration unit (16756) specially developed by Sartorius to the inlet suction port of the AirPort MD8. The individual steps for proper connection of the calibration unit to the AirPort MD8, are described in Section 5 of these Operating Instructions.

The actual calibration can begin once the calibration unit has been properly connected to the AirPort MD8 (Fig. 10) and the service menu on AirPort MD8 has been activated.

Fig. 10

7.2 Activating the Service and Calibration Menu on the AirPort MD8

Note on Activating the Service Menu

To prevent improper operation, and to protect the factory-set default calibration data, a specific key combination was selected for activating the service menu to make it more difficult for the user to access.

In principle, changes to the pre-set calibration should only be performed by trained Sartorius service technicians or if the user is in possession of a suitable calibration unit, e.g. the Sartorius calibration unit 16756.

a) Activating the Service Menu

You can call up the service menu as follows: Press and hold down the $\downarrow\uparrow$ buttons. Now press the I/O button within 1 second, and then let them go.

The display reads:

Service Program

b) Activating the Calibration Menu Item

Go to the calibration menu item by accessing the setting

Service Program

by pressing the $\downarrow\uparrow$ keys 4 times. The following settings will appear in the display, e.g..

Calibration
F-Actual=0 K=089

F-Actual = means the actual value of the current frequency at which the impeller wheel inside the AirPort MD8 scans a fork light barrier.

K-value = the K-value (correction factor, dimensionless) is used to change the air flow rate.

7.3 Calibration Procedure

After connecting the calibration unit and pressing the Start button on the calibration unit, the message "Sartorius MD8 Calibrating System, start air sampler" appears in the display.

On the AirPort MD8, start the actual calibration now by accessing the menu item

Calibration
F-Actual=0 K=089

and pressing the On|Start|Stop button (I/O button):

Press the On|Start|Stop-button. AirPort MD8 regulates the air flow to the previously set flow rate (30, 40 or 50 l/min).

After adjusting the flow rate to either 30, 40 or 50 l/min, now compare actual value indicated on the calibration unit (e.g. 3.0 m³/h or 50 l/min) with the selected nominal value 50 l/min on the AirPort MD8.

If a deviation from the nominal value is noted, you must change the K-value by pressing the \downarrow or \uparrow key until the actual value indicated on the calibration unit agrees with the nominal value.

If the K-value is higher, that means that the air flow rate has increased, a drop in the K-value indicates a reduction in the flow rate.

After this the calibration is completed.

You can save the values after calibration, by pressing the $\downarrow\uparrow$ keys to access the "Service Program" display and pressing On|Start|Stop button.

Note:

The impeller wheel anemometer works within the measuring range of 20 to 80 l/min linear, i.e. only has to be calibrated to one of the 3 air flow rates (30, 40, 50 l/min).

Saving and ending the Service Program is possible from every position.

To save or end, proceed as follows: Use the $\downarrow\uparrow$ keys to access the "Service Program" Display and press On|Start|Stop.

Evaluation:

The flow rate deviation should not exceed 6.66% of the max. value (50 l/min), i.e. ± 3.3 l/min. If the deviation after calibration is greater, it could be caused by an error or a malfunction of the device. In this case, inform the Sartorius Service Department.

Service:

Regular servicing of the calibration unit by the Sartorius service department guarantees constant measuring accuracy.

Sartorius offers service contracts with your choice of maintenance intervals.

Repairs:

Disconnect a defective unit from the power outlet! Repair work may be performed by trained service technicians authorized by Sartorius only. Unauthorized persons may not, under any circumstances, open the air sampler or the calibration unit to re-calibrate it.

Any attempt by untrained persons to perform repairs may lead to hazards for the user.

Always have defective or damaged cables or cable connections (power supply) replaced as an assembled unit.

Cleaning:

Observe the IP protection rating! Make sure that no dust or liquid enters the calibration unit. Prior to any cleaning, servicing or repair work, unplug the unit's power cord from the wall outlet.

Cleaning stainless steel surfaces:
The stainless steel housing can be cleaned at regular intervals. Clean any stainless steel parts with a moistened towel. Commercial cleaning agents suitable for stainless steel can be used without any hazard. For disinfecting surfaces, use commercial disinfectants, suitable for stainless steel.

Do not use any detergents on the stainless steel housing that contain sodium hydroxide solution, acetic acid, sulfuric acid or citric acid. The use of steel wool pads for cleaning is prohibited.

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials

that can be used as secondary raw materials.

The equipment, including accessories and empty non-rechargeable and rechargeable batteries, does not belong in your regular household waste; such equipment is manufactured from high-grade materials and can be recycled and reused. The European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) requires that electrical and electronic equipment be collected and disposed of separately from other unsorted municipal waste with the aim of recycling it. The crossed out waste bin symbol indicates that separate collection is required.

In Germany and many other countries, Sartorius AG takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. Such equipment may not be thrown out with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators. Please contact our on-site service technicians or our Service Center in Goettingen:

Sartorius AG

Service Center
Weender Landstrasse 94-108
37075 Goettingen, Germany

Tel.: +49.551.308.3333 | 4444
Fax: +49.551.308.3730

E-Mail:

Jan.Naumann@sartorius.com

For equipment disposal in the following countries of the European Economic Area (EEA), please contact your local Sartorius subsidiary or distributor:

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Sartorius Stedim Systems GmbH
Schwarzenberger Weg 73-79
34212 Melsungen
Phone +49.5661.71.3400
Fax +49.5661.71.3702
www.sartorius-stedim.com

France

Sartorius Stedim Biotech S.A.
ZI Les Paluds Avenue de Jouques –
BP 105113781 Aubagne Cedex
Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI Les Paluds Avenue de Jouques –
CS 7105813781 Aubagne Cedex
Phone +33.442.845600
Fax +33.442.846545

Austria

Sartorius Stedim Austria GmbH
Franzosengraben 12
A-1030 Vienna
Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Stedim Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde
Phone +32.2.756.06.80
Fax +32.2.756.06.81

Denmark

Sartorius Stedim Nordic A/S
Hoerskaetten 6D, 1.
DK-2630 Taastrup
Phone +45.7023.4400
Fax +45.4630.4030

Hungary

Sartorius Stedim Hungária Kft
Kagyló u. 5
2092 Budakeszi
Phone +36.23.457.227
Fax +36.23.457.147

Italy

Sartorius Stedim Italy S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)
Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Stedim Netherlands B.V.
Edisonbaan 24
3439 MN Nieuwegein
Phone +31.30.6025080
Fax +31.30.6025099

Poland

Sartorius Stedim Poland Sp. z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn
Phone +48.61.647.38.40
Fax +48.61.879.25.04

Spain

Sartorius Stedim Spain SA
C/Isabel Colbrand 10,
Oficina 70
Polígono Industrial de Fuencarral
28050 Madrid
Phone +34.90.2110935
Fax +34.91.3589623

Switzerland

Sartorius Stedim Switzerland AG
Ringstr. 24 a
8317 Tagelswangen
Phone +41.52.354.36.36
Fax +41.52.354.36.46

U.K.

Sartorius Stedim UK Limited
Longmead Business Park
Blenheim Road, Epsom
Surrey KT19 9 QQ
Phone +44.1372.737159
Fax +44.1372.726171

In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.

Sartorius AG, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination). Before returning any equipment to us, the sender must fill out the attached Non-hazardous Declaration completely and return it to us in advance. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Prior to disposal or scrapping of the equipment, remove any rechargeable or non-rechargeable batteries and, if they are empty, dispose of them in local collection boxes. Only for customers in Germany: On request, Sartorius will have the Foundation of Manufacturers' Mutual Return System for Batteries (GRS, Stiftung Gemeinsames Rücknahmesystem Batterien) provide customers in Germany with battery collection boxes.

Note: In Germany, Sartorius AG is entered in the EAR Registry of the Foundation for Used Electrical and Electronic Equipment (Stiftung Elektro-Altgeräte Register). The WEEE registration number for Germany is DE 49923090.

Non-hazard Declaration

Declaration of Safety and Acceptability as Non-Hazardous to Health

The safety and health of our employees, the legal regulations governing the handling of materials that pose hazards to health and the environment and the regulations governing occupational safety mandate and require that this Non-hazardous Declaration be filled out completely for all equipment returned to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Fax or mail a copy of a completed Non-hazardous Declaration form in advance to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country (for the fax number or address, please refer to our list of countries). We require this information before the equipment | component arrives. An additional copy must accompany the equipment | component. If necessary, the shipping company must be notified.

Incomplete entries or non-compliance with this procedure will automatically lead to considerable delays in processing.

In the event of return, we do not differentiate between

- still utilizable and | or repairable used electrical and electronic equipment and no longer utilizable and | or repairable used electrical and electronic equipment or
- between waste that is subject to monitoring by the supervisory authorities and waste that is not subject to any special monitoring by the supervisory authorities (used electrical and electronic equipment).

Model: _____ Serial No.: _____

I | we guarantee that the following measures have been taken:

- ☐ The equipment | component has been freed of hazardous materials to ensure that persons involved with the handling | repair are not exposed to any hazard or danger whatsoever.
- ☐ The equipment has been packaged for safety and the complete labeling has been affixed to it.
- ☐ The shipping agent has been notified (as prescribed) about the hazardous nature of the shipment.

The person sending in the equipment | component is hereby notified that said person shall be liable to Sartorius, its affiliates, subsidiaries, dealers and distributors as well as to any third parties – in particular to any of employees of Sartorius, its affiliates, its subsidiaries, dealers or distributors involved with the handling | repair of the equipment | component – for any damage, caused by entries that are incomplete or incorrect on account of negligence, gross negligence or willful intent.

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Sartorius Stedim Biotech GmbH
August-Spindler-Str. 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.32 89
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part of
this publication may be reprinted
or translated in any form or by
any means without the prior
written permission of Sartorius
Stedim Biotech GmbH.
The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology, fea-
tures, specifications and design
of the equipment without notice.

Status:
November 2011,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany